



ACCESSORIES

Catheter Repair Clamp

CardioMed Supplies Inc.



CardioMed is pleased to introduce a fast and simple solution for replacing defective pinch clamps on tunneled catheters.

Our catheter repair clamp has uniquely designed grooves that allow for immediate mounting of the repair clamp onto catheter extension tubes. The grooves ensure that the clamp will remain in the same linear position as the original clamp.

- ✓ **Unique grooves allow the repair clamp to be mounted easily and securely**
- ✓ **Guide pin keeps the clamp centered ensuring that the extension tube is pinched evenly**

ORDERING INFORMATION			
ITEM NUMBER	DESCRIPTION	SIZE	PACKAGING
CM-604	REPAIR CLAMP	fits .118 inch - .177 inch OD Tubing	BOX: Six (6) Clamps per Box NON-STERILE
LATEX FREE PRODUCT			
CardioMed Supplies Inc. , 199 St. David Street, Lindsay, ON K9V 5K7 Tel: 705.328.2518 • Toll Free: 800.387.9757 • Fax: 705.328-9747 Web: www.cardiomed.com • Email: mail@cardiomed.com Quality Management System registered to ISO 13485:2003 under CMDCAS, Annex V of MDD 93/42/EEC, FDA			Over 27 years of medical excellence...